


B45 Medical Keyboard Upgrade Guide


1. Click “1026 Bastron_Keyboard_V3.1.4” to install the software:




2. Open the software.




3. Connect the keyboard through the USB mode to the computer.
Let the keyboard be in the state of USB mode connection.




4. Press the cleaning lock button until the indicator lights on.




5. Then the software interface will be changed as below,
Please click “Advance” key




6. In the pop-up dialog, select “Y” Key.




7. Please click “File” key




8. Find and open the “segment1_B45-1-nRF51822-APP_B4501HV02-20180818V23.bin” file from upgrade package.




9. Please click “Update” key



10. When the interface shows 100%.



11. Please click “Done” key



12. Please do not click on the content that is not related to upgrade in the software.
Finally turn off the software, then finish the upgrade.